


OPA

High-level disinfectant for medical devices, heat-sensitive dental & endoscopic instruments

HISEPT OPA is a high-level disinfectant for reusable medical devices, patient-care equipment, heat-sensitive dental and endoscopic instruments. This product contains ortho-phthalaldehyde 0.575% as a disinfection agent, with fast-acting and broad-spectrum antimicrobial activity. It eliminates virtually all recognized pathogenic microorganisms on medical devices, including all steam and heat-sensitive semi-critical instruments.



Characteristics

- Contains 0.575% Ortho-Phthalaldehyde
- Broad-spectrum antimicrobial activity
- Ready to use without the need of any activator
- Reusable for up to 14 days

Product Specification

| | |
|------------|---------------------------------|
| Color/Form | Clear & colorless solution |
| pH | 6 - 8 |
| Scent | Characteristic odor of aldehyde |
| Shelf Life | 24 months |



THIS IS A HISEPT PRODUCT, WITH
**HIGH-LEVEL
SAFETY &
PROTECTION**

Ingredients

- 100 g HISEPT OPA contains:
- 0.575 g ortho-phthalaldehyde

Auxiliary ingredients:

- Chelating agents
- Buffering agent
- Corrosion inhibitors
- Deionized water

Microbiological Effectiveness

Meets European Standards

Bactericidal Activity

According to EN 13727 Staphylococcus. aureus/ ATCC 6538
(contact time: 10 minutes) Enterococcus. hirae/ ATCC 10541
Pseudomonas. aeruginosa/ ATCC 15442
Escherichia. coli/ ATCC 10536

According to EN 14561 Staphylococcus. aureus/ ATCC 6538
(contact time: 10 minutes) Enterococcus. hirae/ ATCC 10541
Pseudomonas. aeruginosa/ ATCC 15442

Yeasticidal Activity

According to EN 13624 Candida. albicans/ ATCC 10231
(contact time: 10 minutes)

According to EN 14562 Candida. albicans/ ATCC 10231
(contact time: 10 minutes)

Fungicidal Activity

According to EN 13624 Aspergillus. niger/ ATCC 16404
(contact time: 10 minutes)

According to EN 14562 Aspergillus. niger/ ATCC 16404
(contact time: 15 minutes)

Sporicidal Activity

According to EN 13704 Bacillus. subtilis/ ATCC 6633
(contact time: 60 minutes)



Directions for use

Immerse clean device completely into the ready-to-use HISEPT OPA, filling all lumens and eliminating bubbles, then cover the tank for at least 10 minutes. Remove device from the solution and rinse thoroughly with clean water.

Compatibility

HISEPT OPA is compatible with most materials such as stainless steel, copper, polyethylene, rubber, PVC, polyurethane and silicon.


2/2 Latest revision: Sep 1, 2023

Packagings

| | |
|-------------------------|----------|
| 3 canisters per carton, | 5 L each |
| 6 bottles per carton, | 1 L each |

Hazard Rating

| Ready-to-use | NFPA |
|--------------|------|
| Health | 2 |
| Flammability | 0 |
| Reactivity | 0 |



Warning

- The disinfection tank shall always be covered with a lid. If the solution is kept in a closed container, it will be reusable for up to 14 days.
- In case of contact with skin and eyes, rinse immediately with plenty of water and seek medical advice.
- Keep out of reach of children and unaware persons.
- Store between + 5 °C and + 35 °C and away from direct sunlight
- Contact with Ortho-Phthalaldehyde solution may stain exposed skin or clothing. It is recommended to use nitrile or butyl rubber gloves. When using latex gloves, the user should use double gloves or change single gloves frequently.



Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC

HISEPT[®]
HIGH-LEVEL SAFETY & PROTECTION

Middle East Office
Behban Pharmed Lotus Co.
Address : Unit 113, 12 Floor, Bldg. B, Kaveh
Complex, Azadi Ave., Tehran, Iran
Tel: (+9821) 66581294-5/66581269

info@hisept.com

World Locations: Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad