



RAPID

Fast- acting disinfectant for medical equipment and surfaces

HISEPT RAPID is a ready-to-use and rapid disinfectant for medical equipment and surfaces. Ultra-fast-acting formulation besides broad-spectrum bactericidal and fungicidal properties, combined with its excellent stability and environmentally friendly characteristics, have made it the best choice for disinfection of all alcohol-resistant medical equipment while being an ideal rapid surface disinfectant.



Characteristics

- Ready-to-use disinfectant
- Fast-acting formulation to kill germs in 3 minutes
- Contains no dangerous phenols, aldehydes or chlorine
- User-friendly design of sprayer head

Product Specification

Color/Form	Colorless to pale yellow clear solution
pH	7-10
Scent	Alcohol-like odor
Shelf Life	24 months



THIS IS A HISEPT PRODUCT, WITH
**HIGH-LEVEL
SAFETY &
PROTECTION**

Ingredients

100 g HISEPT RAPID contains:

- 0.6 g N-(3-aminopropyl)-N-dodecylpropane-1, 3-diamine
- 65 g Ethanol

Auxiliary ingredients:

- Surfactant
- pH-Adjusting agent
- Essence
- Deionized water

Microbiological Effectiveness

Meets European Standards

Bactericidal Activity

According to EN 1276 (contact time: 3 minutes)	<i>Pseudomonas. aeruginosa</i> / ATCC 15442 <i>Escherichia. coli</i> / ATCC 10536 <i>Staphylococcus. aureus</i> / ATCC 6538 <i>Enterococcus. hirae</i> / ATCC 10541
According to EN 13727 (contact time: 3 minutes)	<i>Staphylococcus. aureus</i> / ATCC 6538 <i>Enterococcus. hirae</i> / ATCC 10541 <i>Pseudomonas. aeruginosa</i> / ATCC 15442
According to EN 13697 (contact time: 3 minutes)	<i>Staphylococcus. aureus</i> / ATCC 6538 <i>Enterococcus. hirae</i> / ATCC 10541 <i>Pseudomonas. aeruginosa</i> / ATCC 15442 <i>Escherichia. coli</i> / ATCC 10536

Mycobactericidal Activity

According to EN 14348 (contact time: 3 minutes)	<i>Mycobacterium. avium</i> / ATCC 15769 <i>Mycobacterium. terrae</i> / ATCC 15755
---	---

Yeasticidal Activity

According to EN 13624 (contact time: 3 minutes)	<i>Candida. albicans</i> / ATCC 10231
According to EN 13697 (contact time: 3 minutes)	<i>Candida. albicans</i> / ATCC 10231

Fungicidal Activity

According to EN 13624 (contact time: 3 minutes)	<i>Aspergillus. niger</i> / ATCC 16404
According to EN 13697 (contact time: 3 minutes)	<i>Aspergillus. niger</i> / ATCC 16404

Directions for use

Spray directly to the surface or device or spray to a clean linen cloth. Wipe and wait for at least 3 minutes. Do not spray directly to electronic devices.

Compatibility

HISEPT RAPID is compatible with most materials such as stainless steel, copper, plastic, rubber and silicon.

2/2 Latest revision: May 1, 2022

Packagings

6 bottles per carton,	1 L each
20 bottles per carton,	500 mL each
30 bottles per carton,	250 mL each
48 bottles per carton,	125 mL each

Hazard Rating

Ready-to-use	NFPA
Health	1
Flammability	3
Reactivity	1



Warning

- Keep out of reach of children and unaware persons.
- Highly flammable. Keep away from heat, spark and high-voltage electric current.
- Store between + 5 °C and + 35 °C and away from direct sunlight.



Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC



HISEPT[®]
HIGH-LEVEL SAFETY & PROTECTION

Middle East Office
Behban Pharmed Lotus Co.
Address : Unit 113, 12 Floor, Bldg. B, Kaveh
Complex, Azadi Ave., Tehran, Iran
Tel: (+9821) 66581294-5/66581269

info@hisept.com

World Locations: Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad