



FOAM

Ready- to- use, cleaner and disinfectant for medical equipment and surfaces

HISEPT FOAM is a ready-to-use and fast-acting foam disinfectant designed for alcohol-sensitive equipment and surfaces. Alkylamine is a powerful antimicrobial agent with potent cidal action on a wide variety of gram-negative and gram-positive bacteria as well as on viruses. HISEPT Foam is cost-effective and suitable for hard-to-reach surfaces and objects, causes no corrosion, provides potent lethal activity against microbial population which prevents resistance.



Characteristics

- Ready-to-use foam cleaner and disinfectant
- Contains no dangerous phenols, aldehydes or chlorine
- Cost-effective
- Excellent cleaning performance

Product Specification

Color/Form	Colorless to yellow clear solution
pH	7-10
Scent	Pleasant odor
Shelf Life	24 months



THIS IS A HISEPT PRODUCT, WITH
**HIGH-LEVEL
SAFETY &
PROTECTION**

Ingredients

100 g HISEPT FOAM contains:

- 0.37 g N-(3-aminopropyl)-N-dodecylpropane-1, 3-diamine

Auxiliary ingredients:

- Ethoxylated alcohol
- pH-Adjusting agent
- Essence
- Deionized water

Microbiological Effectiveness

Meets European Standards

Bactericidal Activity

According to EN 13727 (contact time: 5 minutes) Staphylococcus. aureus/ ATCC 6538
Enterococcus. hirae/ ATCC 10541
Pseudomonas. aeruginosa/ ATCC 15442

According to EN 13697 (contact time: 5 minutes) Staphylococcus. aureus/ ATCC 6538
Enterococcus. hirae/ ATCC 10541
Pseudomonas. aeruginosa/ ATCC 15442
Escherichia. coli/ ATCC 10536

Yeasticidal Activity

According to EN 13624 (contact time: 5 minutes) Candida. albicans/ATCC 10231

According to EN 13697 (contact time: 5 minutes) Candida. albicans/ATCC 10231

Virucidal Activity

According to EN 14476 (contact time: 3 minutes) Adenovirus type 5
Herpes simplex virus (HSV-1)

Packagings

20 bottles per carton,

500 mL each

30 bottles per carton,


250 mL each

48 bottles per carton,

125 mL each

Hazard Rating

Ready-to-use	NFPA
Health	1
Flammability	0
Reactivity	0



Warning

- Keep out of reach of children and unaware persons.
- Store between +5 °C and +35 °C and away from direct sunlight.



Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC



Directions for use

Spray directly to the surface or device or spray to a clean linen cloth. Wipe and wait for at least 3 minutes. Do not spray directly to electronic devices.

Compatibility

HISEPT FOAM is compatible with most materials such as stainless steel, copper, plastic, rubber and silicon.

2/2 Latest revision: April 1, 2023

HISEPT[®]
HIGH-LEVEL SAFETY & PROTECTION

Middle East Office
Behban Pharmed Lotus Co.
Address : Unit 113, 12 Floor, Bldg. B, Kaveh
Complex, Azadi Ave., Tehran, Iran
Tel: (+9821) 66581294-5/66581269

info@hisept.com

World Locations: Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad