

MEDICLEAN 116

Acidic neutralizer for automated reprocessing of surgical instruments

HISEPT MEDICLEAN 116 is an acidic detergent for the neutralization of alkaline residues on surgical instruments. This product can strongly remove limescale, rust stains and mineral deposits of the instruments.

Characteristics

- Strong remover of alkaline and mineral residues
- Highly concentrated and cost-effective
- Restores stainless steel surface finish

Product Specification

Color/Form	Clear & colorless solution
pH (0.5 %)	2 - 4
Scent	Odorless
Shelf Life	24 months



THIS IS A HISEPT PRODUCT, WITH
**HIGH-LEVEL
SAFETY &
PROTECTION**

Ingredients

HISEPT MEDICLEAN 116 contains:

- Organic acids
- Deionized water

Mirobiological Specification

According to	
USP 40	Bacterial count <10 CFU/10g
	Yeast & mold total count < 10 CFU/10g
	No growth of Staphylococcus. aureus
	No growth of Pseudomonas. aeruginosa
	No growth of Escherichia. coli
	No growth of Fungi



Directions for use

This product should be diluted based on the washing machine instructions. Please refer to the manufacturer's recommendations.

Compatibility

HISEPT MEDICLEAN 116 is compatible with most materials such as stainless steel, copper, plastic, rubber and silicon.



2/2 Latest revision: May 1, 2022

Packagings

2 canisters per carton,	10 L each
3 canisters per carton,	5 L each

Hazard Rating

Concentrate	NFPA
Health	1
Flammability	0
Reactivity	0



Warning

- Apply the safety accessories such as work gloves, boots, mask, and work clothes during using & dilution.
- If it is in contact with skin and eyes, rinse immediately with plenty of water.
- Keep in the original container and sealed .
- Keep out of the combustibile materials and alkaline compounds.
- Keep out of reach of children and keep in a ventilated area.



Certifications

- >> ISO 9001: 2015
- >> ISO 13485: 2016
- >> GMP
- >> FSC



HISEPT[®]
HIGH-LEVEL SAFETY & PROTECTION

Middle East Office
Behban Pharmed Lotus Co.
Address : Unit 113, 12 Floor, Bldg. B, Kaveh
Complex, Azadi Ave., Tehran, Iran
Tel: (+9821) 66581294-5/66581269

info@hisept.com

World Locations: Düsseldorf, Kuala Lumpur, Istanbul, Tehran, Baghdad